

Convoid Rodeo Tilt w/out Seating - HCPCS E1234

Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____

Contact: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

Rodeo Tilt Standard Features:

- 5° to 45° of Adjustable Tilt
- Silver Powder Coated Frame
- Adjustable One-Piece Push Handle
- 2-Point Positioning Belt
- Elevating Legrest / One-Piece, Angle Adjustable Footplate
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tires
- Quick Release Rolko Rear Wheels
- Single Foot Activated Wheel Lock Bar
- Seat Depth Growth of 5"
- Transit Anchor Brackets
- Made in the USA

Standard Models:

	Price
<input type="checkbox"/> Rodeo Transit w/out Seating RD10TWOS - 904967	\$2,603
<input type="checkbox"/> Rodeo Transit w/out Seating RD12TWOS - 904968	\$2,603
<input type="checkbox"/> Rodeo Transit w/out Seating RD14TWOS - 904969	\$2,603
<input type="checkbox"/> Rodeo Transit w/out Seating RD16TWOS - 904970	\$2,603

	Price	HCPCS	RD10	RD12	RD14	RD16
Seating Module						
<input type="checkbox"/> Depth Adjustable Solid Seat Pan (Manufacturer Installation Required)	\$150	E2293	904971	904972	904973	904974
Transit Required Items						
<input type="checkbox"/> Transit Option (Manufacturer Installation Required)	\$0		904341B	903167B	903168B	903169B
<input type="checkbox"/> Using other manufacturer's Head Support	\$0		995681	995681	995681	995681
<input type="checkbox"/> 2-Point Positioning Belt	\$0	E0978	903931B	903931B	904559B	904559B
Transit Options						
<input type="checkbox"/> Q'Straint Transit Lap Belt	\$90	E0978	903094	903094	903094	903094
<input type="checkbox"/> Sure-Lok Transit Lap Belt	\$90	E0978	903232	903232	903232	903232
<input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Small ^	\$151	E0950	903669	903669		
<input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Medium ^	\$151	E0950			903670	903670
^ Crash tested Foam Tray with vinyl cover. Flip-Up Armrests required.						
Foot & Leg Positioning						
<input type="checkbox"/> Individual Angle Adjustable Footplates	\$179	K0040		904957	904958	904959
<input type="checkbox"/> Full Padded Footbox (One-Piece Footplate Required)	\$325	K0108	904953	904954	904955	904956
<input type="checkbox"/> Legrest Extension - 2" *	\$55		904095	904095	904096	904097
* 2" Legrest Extension interferes with caster swivel when seat is in upright position. Seat must be tilted back at least 5 degrees OR Legrests must be elevated to allow caster swivel.						
Frame Option						
<input type="checkbox"/> Recline Option (Manufacturer Installation Required)	\$461	E1225	904342	904343	904344	904345
<input type="checkbox"/> Height Adjustable Flip-Up Armrests (pair)	\$228	E0973	904346	904346	904346	904346
<input type="checkbox"/> Solid Side Panels	\$211	K0108	904388	902672	902652	903145
<input type="checkbox"/> Mesh Side Panels	\$148	K0108	904389	902680	902681	903146
Tire Options						
<input type="checkbox"/> 12.5" Rear Solid Knobby Tire	\$199		904117	904117	904117	904117
<input type="checkbox"/> 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tire	\$299		904118	904118	904118	904118
Accessories						
<input type="checkbox"/> Headrest Cover (Canopy)	\$228		904366	902679	902659	903152
<input type="checkbox"/> Extended Headrest Cover - No Windows (Canopy)	\$250		904728	904729	904730	904731
<input type="checkbox"/> Extended Headrest Cover - with Windows (Canopy)	\$295		904744	904745	904746	904747
<input type="checkbox"/> Upper Extremity Support Surface (Clear Tray) ^	\$111	E0950	904367	903405	903401	903402
<input type="checkbox"/> Adj. Removable Hardware (Required for Clear Tray)	\$227	E1028	903369	903369	903369	903369
<input type="checkbox"/> Medical Necessity Storage Basket	\$295		904476	904476	904477	904477
<input type="checkbox"/> Under Seat Storage Basket - Tilt-in-Space	\$75		904529	904529	904530	904530
<input type="checkbox"/> IV Pole (height adjustable and collapsible) ²	\$150	K0105	904566	904566	904566	904566
<input type="checkbox"/> Oxygen Tank Bag	\$179	E2208	902870	902870	902870	902870
<input type="checkbox"/> Utility Bag	\$38		904531	904531	904531	904531
<input type="checkbox"/> Attendant Hand Brakes (pair) ¹	\$236				904370	904371
<input type="checkbox"/> Curb Tipper	\$98		905310	905310	905310	905310
<input type="checkbox"/> Caster Lock for 2" Wide Front Tires ¹	\$94	K0073	904650	904650	904650	904650
<input type="checkbox"/> Rear Anti-Tip Tubes (Pair)*	\$115	E0971	904649	904649	904649	904649
<input type="checkbox"/> Rain Cover Size 1 (Canopy Req)	\$150		904939	904939	904939	
<input type="checkbox"/> Rain Cover Size 2 (Canopy Req) (RD14 Extended Canopy Req)	\$150				904940	904940
<input type="checkbox"/> Rain Cover Size 3 (Extended Canopy Req)	\$150					905283
<input type="checkbox"/> Mosquito Net Size 1 (Canopy req)	\$150					
<input type="checkbox"/> Mosquito Net Size 2 (Canopy req) (RD14 Extended Canopy Req)	\$150		904655	904655	904655	
<input type="checkbox"/> Mosquito Net Size 3 (Extended Canopy req)	\$150				904656	904656
<input type="checkbox"/> Flip-Up Armrests and Hardware required.						905319

^ Flip-Up Armrests and Hardware required.
 *Highly recommended for users weighing 50lbs or less.
¹ Manufacturer Installation Required.
² Anti-tip Tubes Required

The Rodeo Transit option includes four manufacturer installed red wheelchair transportation anchors. The Rodeo Transit model has been crash tested and conforms to RESNA WC19 standards. Anthropomorphic Test Device weights are RD10T, RD12T 66lbs/30kg, RD14T 100lbs/45kg, and RD16T 140lbs/64kg.

The Convaid Rodeo Mobility Base was tested with the Rodeo Seating Module, equipped with a head support and wheelchair-anchored 2-point pelvic lap belt. Do not use the foot positioners on the Convaid Rodeo during transit in a motor vehicle. The performance and safety of accessories and/or other seating options not tested by Convaid is unknown, and no warranty or claim about their performance is or can being made.

It is Convaid's position that the wheelchair user should be transferred from the wheelchair into appropriate vehicle seating whenever possible. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair, therefore Convaid cannot recommend a transit system.

The WC19 standard is specifically designed for the wheelchair to be secured in a forward facing configuration and not in any other direction. For more information on RESNA WC-19 and WC-20 standards please visit <http://www.ercwts.pitt.edu/WC19.html>.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245

2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166

Prices in USD & subject to change, ex works Torrance, California.

Effective: 8/24/2020

www.convaid.com

OF0110-ENG-REV.02