

# Convoid Rodeo Tilt w/out Seating - HCPCS E1234

## Retail Price List / Order Form

Account #: \_\_\_\_\_

Purchase Order #: \_\_\_\_\_

Date: \_\_\_\_\_

Contact: \_\_\_\_\_

Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

### Rodeo Tilt Standard Features:

- 5° to 45° of Adjustable Tilt
- Silver Powder Coated Frame
- Adjustable One-Piece Push Handle
- Elevating Legrest / One-Piece, Angle Adjustable Footplate
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tires
- Quick Release Rolko Rear Wheels
- Single Foot Activated Wheel Lock Bar
- Seat Depth Growth of 5"
- Made in the USA

#### Standard Models:

	<u>Price</u>
<input type="checkbox"/> Rodeo w/out Seating RD10WOS - 904963	<b>\$2,318</b>
<input type="checkbox"/> Rodeo w/out Seating RD12WOS - 904964	<b>\$2,318</b>
<input type="checkbox"/> Rodeo w/out Seating RD14WOS - 904965	<b>\$2,318</b>
<input type="checkbox"/> Rodeo w/out Seating RD16WOS - 904966	<b>\$2,318</b>

<u>Seating Module</u>	<u>Price</u>	<u>HCPCS</u>	<u>RD10</u>	<u>RD12</u>	<u>RD14</u>	<u>RD16</u>
<input type="checkbox"/> Depth Adjustable Solid Seat Pan (Manufacturer Installation Required)	<b>\$150</b>	<b>E2293</b>	904971	904972	904973	904974

<u>Transit Required Items</u>	<u>Price</u>	<u>HCPCS</u>	<u>RD10</u>	<u>RD12</u>	<u>RD14</u>	<u>RD16</u>
<input type="checkbox"/> Transit Option (Manufacturer Installation Required)	<b>\$200</b>		904341	903167	903168	903169
<input type="checkbox"/> Using other manufacturer's Head Support	<b>\$0</b>		995681	995681	995681	995681
<input type="checkbox"/> 2-Point Positioning Belt	<b>\$45</b>	<b>E0978</b>	903931	903931	904559	904559

<u>Transit Options</u>	<u>Price</u>	<u>HCPCS</u>	<u>RD10</u>	<u>RD12</u>	<u>RD14</u>	<u>RD16</u>
<input type="checkbox"/> Q'Straint Transit Lap Belt	<b>\$90</b>	<b>E0978</b>	903094	903094	903094	903094
<input type="checkbox"/> Sure-Lok Transit Lap Belt	<b>\$90</b>	<b>E0978</b>	903232	903232	903232	903232
<input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Small ^	<b>\$151</b>	<b>E0950</b>	903669	903669		
<input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Medium ^	<b>\$151</b>	<b>E0950</b>			903670	903670

<u>Foot &amp; Leg Positioning</u>	<u>Price</u>	<u>HCPCS</u>	<u>RD10</u>	<u>RD12</u>	<u>RD14</u>	<u>RD16</u>
<input type="checkbox"/> Individual Angle Adjustable Footplates	<b>\$179</b>	<b>K0040</b>		904957	904958	904959
<input type="checkbox"/> Full Padded Footbox (One-Piece Footplate Required)	<b>\$325</b>	<b>K0108</b>	904953	904954	904955	904956
<input type="checkbox"/> Legrest Extension - 2" *	<b>\$55</b>		904095	904095	904096	904097
<input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair)	<b>\$156</b>	<b>K0108</b>	903415	903415	903415	903415
<input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair)	<b>\$156</b>	<b>K0108</b>	903416	903416	903416	903416
<input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair)	<b>\$156</b>	<b>K0108</b>	903417	903417	903417	903417

\* 2" Legrest Extension interferes with caster swivel when seat is in upright position.  
Seat must be tilted back at least 5 degrees OR Legrests must be elevated to allow caster swivel.

<u>Frame Option</u>	<u>Price</u>	<u>HCPCS</u>	<u>RD10</u>	<u>RD12</u>	<u>RD14</u>	<u>RD16</u>
<input type="checkbox"/> Recline Option (Manufacturer Installation Required)	<b>\$461</b>	<b>E1225</b>	904342	904343	904344	904345
<input type="checkbox"/> Height Adjustable Flip-Up Armrests (pair)	<b>\$228</b>	<b>E0973</b>	904346	904346	904346	904346
<input type="checkbox"/> Solid Side Panels	<b>\$211</b>	<b>K0108</b>	904388	902672	902652	903145
<input type="checkbox"/> Mesh Side Panels	<b>\$148</b>	<b>K0108</b>	904389	902680	902681	903146

<u>Tire Options</u>	<u>Price</u>	<u>HCPCS</u>	<u>RD10</u>	<u>RD12</u>	<u>RD14</u>	<u>RD16</u>
<input type="checkbox"/> 12.5" Rear Solid Knobby Tire	<b>\$199</b>		904117	904117	904117	904117
<input type="checkbox"/> 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tire	<b>\$299</b>		904118	904118	904118	904118

<u>Canopies</u>	<u>Price</u>	<u>HCPCS</u>	<u>RD10</u>	<u>RD12</u>	<u>RD14</u>	<u>RD16</u>
<input type="checkbox"/> Headrest Cover (Canopy)	<b>\$228</b>		904366	902679	902659	903152
<input type="checkbox"/> Extended Headrest Cover - No Windows (Canopy)	<b>\$250</b>		904728	904729	904730	904731
<input type="checkbox"/> Extended Headrest Cover - with Windows (Canopy)	<b>\$295</b>		904744	904745	904746	904747
Canopy Upholstery (Please Choose One Below)						
<input type="checkbox"/> Royal Blue	<b>\$0</b>		903850	903850	903850	903850
<input type="checkbox"/> Forest Green	<b>\$0</b>		903857	903857	903857	903857
<input type="checkbox"/> Cherry Red	<b>\$0</b>		903855	903855	903855	903855
<input type="checkbox"/> Sassy Purple	<b>\$0</b>		903851	903851	903851	903851
<input type="checkbox"/> Princess Pink	<b>\$0</b>		903759	903759	903759	903759
<input type="checkbox"/> Panther Black	<b>\$0</b>		903854	903854	903854	903854
<input type="checkbox"/> Chocolate Brown	<b>\$0</b>		904189	904189	904189	904189
<input type="checkbox"/> Aqua Blue	<b>\$0</b>		904753	904753	904753	904753
<input type="checkbox"/> Apple Green	<b>\$0</b>		904754	904754	904754	904754

<b>Accessories</b>	<b>Price</b>	<b>HCPCS</b>
<input type="checkbox"/> Upper Extremity Support Surface (Clear Tray) ^	<b>\$111</b>	<b>E0950</b>
<input type="checkbox"/> Adj. Removable Hardware (Required for Clear Tray)	<b>\$227</b>	<b>E1028</b>
<input type="checkbox"/> Medical Necessity Storage Basket	<b>\$295</b>	
<input type="checkbox"/> Under Seat Storage Basket - Tilt-in-Space	<b>\$75</b>	
<input type="checkbox"/> IV Pole (height adjustable and collapsible) <sup>2</sup>	<b>\$150</b>	<b>K0105</b>
<input type="checkbox"/> Oxygen Tank Bag	<b>\$179</b>	<b>E2208</b>
<input type="checkbox"/> Utility Bag	<b>\$38</b>	
<input type="checkbox"/> Attendant Hand Brakes (pair) <sup>1</sup>	<b>\$236</b>	
<input type="checkbox"/> Curb Tipper	<b>\$98</b>	
<input type="checkbox"/> Caster Lock for 2" Wide Front Tires <sup>1</sup>	<b>\$94</b>	<b>K0073</b>
<input type="checkbox"/> Rear Anti-Tip Tubes (Pair)*	<b>\$115</b>	<b>E0971</b>
<input type="checkbox"/> Rain Cover Size 1 (Canopy Req)	<b>\$150</b>	
<input type="checkbox"/> Rain Cover Size 2 (Canopy Req) (RD14 Extended Canopy Req)	<b>\$150</b>	
<input type="checkbox"/> Rain Cover Size 3 (Extended Canopy Req)	<b>\$150</b>	
<input type="checkbox"/> Mosquito Net Size 1 (Canopy req)	<b>\$150</b>	
<input type="checkbox"/> Mosquito Net Size 2 (Canopy req) (RD14 Extended Canopy Req)	<b>\$150</b>	
<input type="checkbox"/> Mosquito Net Size 3 (Extended Canopy req)	<b>\$150</b>	

RD10	RD12	RD14	RD16
904367	903405	903401	903402
903369	903369	903369	903369
904476	904476	904477	904477
904529	904529	904530	904530
904566	904566	904566	904566
902870	902870	902870	902870
904531	904531	904531	904531
		904370	904371
905310	905310	905310	905310
904650	904650	904650	904650
904649	904649	904649	904649
904939	904939	904939	
		904940	904940
			905283
904655	904655	904655	
		904656	904656
			905319

^ *Flip-Up Armrests and Hardware required.*  
 \**Highly recommended for users weighing 50lbs or less.*  
<sup>1</sup> *Manufacturer Installation Required.*  
<sup>2</sup> *Anti-tip Tubes Required.*

The Rodeo Transit option includes four manufacturer installed red wheelchair transportation anchors. The Rodeo Transit model has been crash tested and conforms to RESNA WC19 standards. Anthropomorphic Test Device weights are RD10T, RD12T 66lbs/30kg, RD14T 100lbs/45kg, and RD16T 140lbs/64kg.

The Convoid Rodeo Mobility Base was tested with the Rodeo Seating Module, equipped with a head support and wheelchair-anchored 2-point pelvic lap belt. Do not use the foot positioners on the Convoid Rodeo during transit in a motor vehicle. The performance and safety of accessories and/or other seating options not tested by Convoid is unknown, and no warranty or claim about their performance is or can being made.

It is Convoid's position that the wheelchair user should be transferred from the wheelchair into appropriate vehicle seating whenever possible. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair, therefore Convoid cannot recommend a transit system.

The WC19 standard is specifically designed for the wheelchair to be secured in a forward facing configuration and not in any other direction. For more information on RESNA WC-19 and WC-20 standards please visit <http://www.rercwts.pitt.edu/WC19.html>.

*HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.*

**Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245**  
 Prices in USD & subject to change, ex works Torrance, California.

Effective: 11/6/2020

Convoid  
 2830 California Street  
 Torrance, CA 90503  
 1 844-US-MOBILITY (844-876-6245)  
 www.etac.us.com

