

# R82 Quokka™

## USD Retail Price List/ Order Form


Account #: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Date: \_\_\_\_\_  
 Contact: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email	Phone:

**US Configuration Options**


**Effective: 10.20.20 Rev 4**

Quokka Child Restraint System



- A (Back Height) \_\_\_\_\_
- B (Arm Height) \_\_\_\_\_
- C (Lower Leg Length) \_\_\_\_\_
- D (Seat Depth) \_\_\_\_\_
- E (Chest Width) \_\_\_\_\_
- F (Seat Width) \_\_\_\_\_
- G (Height: Seat to Top of Head) \_\_\_\_\_
- H (Shoulder Width) \_\_\_\_\_
- Age \_\_\_\_\_
- Weight \_\_\_\_\_


**Measurements**



**Quokka  
Inches**

Seat Width Maximum (A)	11"
Seat Depth (B)	11"
Height of Shoulder Strap Slots (C)	7 ½", 9 ½", 11 ½", 13 ½"
Back Height (Top of Head) Maximum (D)	21"
Head Support Height Range Minimum	13 ½"
Head Support Height Range Maximum	19 ½"
Head Support Height	6 ½"
Exterior Shell Width (W), Depth (D), Height (H)	16" x 23" x 23 ½"
Interior Shell Width (W), Depth (D), Height (H)	11" x 11" x 21"
Minimum User Weight	5 lbs.
Maximum User Weight	39.6 lbs. (including cast weight)
Maximum User Height	41"
Product Weight	15.8 lbs.

## Quokka Standard Features



- Quokka Seat, One Size, Black Shell
- Head Support (Red)
- 5-Point Safety Harness
- Harness Padded Covers
- Seat Cover (Red/Black Trim)
- Spica Positioning Wedge
- LATCH System (Lower Anchors and Tether Strap)

Description	Size	Part No.	USD Price
<input type="checkbox"/> Quokka Seat, (Red/Black Trim)	One Size	905371-2	\$700

Accessories	Description	Measurements	Size	Part No.	USD Price
<input type="checkbox"/>	Seat Saver, Incontinence Pad		One Size	907029	\$98

Spare Parts	Description	Measurements	Size	Part No.	USD Price
<input type="checkbox"/>	Harness Buckle & Strap * <i>*Replacement for standard part</i>	6"	Standard	905386	\$20
<input type="checkbox"/>	Harness Buckle & Strap	10"	XL	905387	\$20
<input type="checkbox"/>	Replacement Seat Cover		One Size	905376-2	\$165
<input type="checkbox"/>	Replacement 5-Point Harness		One Size	905378	\$85
<input type="checkbox"/>	Extra Tether Strap		One Size	905375	\$46
<input type="checkbox"/>	Extra Lower Anchor Strap		One Size	905374	\$46
<input type="checkbox"/>	Spica Wedge Positioning System		One Size	905373	\$30
<input type="checkbox"/>	Replacement Headrest		One Size	905377-2	\$50
<input type="checkbox"/>	Replacement Shoulder Protectors 7"			905383	\$43
<input type="checkbox"/>	Replacement Shoulder Protectors 4"			905384	\$43
<input type="checkbox"/>	Belt Locking Clip Set		Pair	905382	\$15

**IMPORTANT: Due to safety regulations and restrictions, no returns or exchanges are accepted of car seat products.**

Prices in USD & is subject to change without notice

**For Quotes and Orders please send completed Order Form to: [sales.us@etac.com](mailto:sales.us@etac.com) or fax: 1-800-760-0537**