

# R82 Wallaroo

## USD Retail Price List/ Order Form


Account #: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Date: \_\_\_\_\_  
 Contact: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email	Phone:

### US Configuration


Effective: 10.30.20 - Rev 7

- Wallaroo Child Restraint System
- Pommel
- Positioning Wedges
- Seat Extension




**A** (Back Height) \_\_\_\_\_  
**B** (Arm Height) \_\_\_\_\_  
**C** (Lower Leg Length) \_\_\_\_\_  
**D** (Seat Depth) \_\_\_\_\_  
**E** (Chest Width) \_\_\_\_\_  
**F** (Seat Width) \_\_\_\_\_  
**G** (Height: Seat to Top of Head) \_\_\_\_\_  
**H** (Shoulder Width) \_\_\_\_\_  
 Age \_\_\_\_\_  
 Weight \_\_\_\_\_

### Measurements



	<b>Wallaroo</b>
	<b>Inches</b>
Seat Width Minimum ( <b>A</b> )	11 <sup>3</sup> / <sub>4</sub> "
Seat Width Maximum ( <b>A1</b> )	13 <sup>3</sup> / <sub>4</sub> "
Seat Depth ( <b>B</b> )	10 <sup>1</sup> / <sub>2</sub> "
Seat Depth Fixed with Extension ( <b>B1</b> )	15"
Height of Shoulder Strap Slots ( <b>C</b> )	11", 13", 15", 17"
Back Height (Top of Head) Maximum ( <b>C1</b> )	26 <sup>3</sup> / <sub>4</sub> "
Head Support Height Range Minimum	11"
Head Support Height Range Maximum	17"
Head Support Height	8"
Exterior Shell Width (W), Depth (D), Height (H)	22" x 22" x 31"
Interior Shell Width (W), Depth (D), Height (H)	15 <sup>1</sup> / <sub>2</sub> " x 14" x 29"
Minimum User Weight	22 lbs.
Maximum User Weight	106 lbs.
Product Weight	27 lbs.

## Wallaroo Standard Features



- Wallaroo Seat, One Size, Black
- Head Support
- 5-Point Positioning Harness with Padded Covers
- Inlay Seat and Headrest Cover (Blue, Red or Grey)
- Buckle Guard
- LATCH System (Lower Anchors and Tether Strap)
- Recline Wedge

Description	Size	Part No.	USD Price
<input type="checkbox"/> Wallaroo Seat, Black with Blue Inlay Cover	One Size	907012-1	\$995
<input type="checkbox"/> Wallaroo Seat, Black with Red Inlay Cover	One Size	907012-2	\$995
<input type="checkbox"/> Wallaroo Seat, Black with Grey Inlay Cover	One Size	907012-4	\$995
<input type="checkbox"/> Wallaroo Seat, Bundle, Black with Blue Inlay Cover**	One Size	907011-1	\$1,345
<input type="checkbox"/> Wallaroo Seat, Bundle, Black with Red Inlay Cover**	One Size	907011-2	\$1,345
<input type="checkbox"/> Wallaroo Seat, Bundle, Black with Grey Inlay Cover**	One Size	907011-4	\$1,345

**\*\*Bundle includes Pommel, Seat Extension, Positioning Wedges, and Support Tray**

Accessories	Description	Measurements	Size	Part No.	USD Price
<input type="checkbox"/>	Pommel, Black (Ea.)		One Size	907013-60	\$50
<input type="checkbox"/>	Seat Extension, Black (Ea.)		One Size	907015-60	\$80
<input type="checkbox"/>	Spica Kit*		One Size	907026-60	\$119
	<b>*Kit includes Spica Insert Cushion, Spica Positioning Wedge, and Spica Cushion Crotch Belt 10"</b>				
<input type="checkbox"/>	Positioning Wedges, Black (Pr.)		One Size	907022-60	\$50
<input type="checkbox"/>	Support Tray, Black (Ea.)		One Size	907014-60	\$170
<input type="checkbox"/>	Seat Saver, Incontinence Pad		One Size	907029	\$98

**IMPORTANT: Due to safety regulations and restrictions, no returns or exchanges are accepted of car seat products.**

Spare Parts	Description	Measurements	Size	Part No.	USD Price
<input type="checkbox"/>	Additional Inlay, Seat & Headrest Cover, Blue (Ea.)		One Size	907018-1	\$125
<input type="checkbox"/>	Additional Inlay, Seat & Headrest Cover, Red (Ea.)		One Size	907018-2	\$125
<input type="checkbox"/>	Additional Inlay, Seat & Headrest Cover, Grey (Ea.)		One Size	907018-4	\$125
<input type="checkbox"/>	Spica Insert Cushion (Ea.)		One Size	907023-60	\$75
<input type="checkbox"/>	Spica Positioning Wedge (Ea.)		One Size	907024-60	\$30
<input type="checkbox"/>	Spica Cushion Crotch Belt 10" (Ea.)		One Size	907025-60	\$20

Prices in USD & is subject to change without notice

**For Quotes and Orders please send completed Order Form to: [sales.us@etac.com](mailto:sales.us@etac.com) or fax: 1-800-760-0537**

Etac  
13137 Bleinheim Lane  
Matthews, North Carolina 28105  
+1 844-US-MOBILITY (844-876-6245)  
[www.etac.us.com](http://www.etac.us.com)

**etac**<sup>®</sup>  
Creating Possibilities