

Convoid Carrot 3 Booster Seat - Child Restraint System

Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____

Contact: _____ Mark For: _____

| | |
|-----------------|-----------------|
| Bill To: | Ship To: |
| Address: | Address: |
| City/State/Zip: | |
| Phone: | City/State/Zip: |
| Email: | Phone: |

Carrot 3 Booster Seat Standard Features:

- Carrot 3 Backrest
- Recline Adjustment (90° - 180°)
- Head Support: 8.5" x 10.5" (H x L)
- 4 Point Positioning Harness
- H-Harness Attaching Bracket
- Booster Seat Base
- Lower Anchors Connector System
- 2" Backrest Extension
- Shoulder Protector Wings

Standard Models:

| | Price |
|--|---------|
| <input type="checkbox"/> Carrot 3 Booster Seat CRT3BLACK - CRT3BS-10 | \$2,750 |
| <input type="checkbox"/> Carrot 3 Booster Seat CRT3BLUE - CRT3BS-20 | \$2,750 |
| <input type="checkbox"/> Carrot 3 Booster Seat CRT3RED - CRT3BS-30 | \$2,750 |

Carrot Booster Seat Recommended Item

- | | Price |
|---|-------|
| <input type="checkbox"/> Tether strap (Black) | \$89 |
| <input type="checkbox"/> Top Securement Strap * | \$53 |
| <input type="checkbox"/> Lower Securement Strap * | \$65 |

*For vehicles not equipped with anchor points.

| Black | Blue | Red |
|-----------|-----------|-----------|
| CT3TS0-10 | CT3TS0-10 | CT3TS0-10 |
| CT3HSF-10 | CT3HSF-10 | CT3HSF-10 |
| CT3LFS-10 | CT3LFS-10 | CT3LFS-10 |

Accessories

- | | Price |
|---|-------|
| <input type="checkbox"/> Booster Support Tray (Crash Tested) | \$290 |
| <input type="checkbox"/> 3D Body Balance Pads L/R Set | \$150 |
| <input type="checkbox"/> Seat Angle Adjustment Wedge (Height: 4.25" / Angle: 20°) | \$113 |
| <input type="checkbox"/> Backrest Extension - 4" (Extends back height to 32.8") | \$91 |
| <input type="checkbox"/> Large Crotch Belt | \$238 |

| Black | Blue | Red |
|-----------|-----------|-----------|
| CRT3BT-10 | CRT3BT-20 | CRT3BT-30 |
| CT3BBP-00 | CT3BBP-00 | CT3BBP-00 |
| CT3SA0-10 | CT3SA0-20 | CT3SA0-30 |
| CT3BE4-10 | CT3BE4-20 | CT3BE4-30 |
| CT3CBX-00 | CT3CBX-00 | CT3CBX-00 |

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245

Prices in USD & subject to change, ex works Torrance, California.


Effective: 12/11/2020

Convoid
2830 California Street
Torrance, CA 90503
1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF2100-ENG-REV.02



CLIENT MEASUREMENTS:



Carrot 3 Base Model Specifications:

| | |
|---------------------------------------|-----------------|
| Seat Width: | 16.5" |
| Seat Depth: | 15" - 16.5" * |
| Backrest Height Adjustment: | 28" - 33" |
| Head Support Height: | 8.5" |
| Backrest Height: | 29" - 34" |
| Backrest Extension Height (Standard): | 2" |
| Backrest Extension Height (Option): | 4" |
| Shoulder Wing Height (Standard): | 6" |
| Shoulder Wing Width Adjustment: | 15.25" - 19.25" |
| Product Weight: | 21 lbs. |
| Recline Adjustment: | 90° - 180° |
| Product Overall Width: | 20" - 22" |

User Measurement Guide:

| | |
|----------------------|----------|
| Minimum User Weight: | 79 lbs. |
| Maximum User Weight: | 165 lbs. |
| Minimum User Height: | 55" |
| Maximum User Height: | 70" |

*Factory Setting

2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166