

# Convoid Rodeo Tilt - HCPCS E1232/E1161

## Retail Price List / Order Form

Account #: \_\_\_\_\_

Purchase Order #: \_\_\_\_\_

Date: \_\_\_\_\_

Contact: \_\_\_\_\_

Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

### Rodeo Tilt Standard Features:

- |   |  |
|---|--|
| <ul style="list-style-type: none"> <li>• 5° to 45° of Adjustable Tilt</li> <li>• Silver Powder Coated Frame</li> <li>• One Piece Height Adjustable Push Handle</li> <li>• One Piece Angle Adjustable Footplate</li> <li>• Seat Depth Growth of 5"</li> <li>• 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire</li> <li>• Quick Release Rolko Rear Tires</li> <li>• Single Foot Activated Wheel Lock Bar</li> <li>• HCPCS code E1161 is only valid on size 16 (RD16)</li> </ul> | <ul style="list-style-type: none"> <li>• Transit Anchor Brackets</li> <li>• H-Harness with Padded Covers</li> <li>• 3-Pt Positioning Belt</li> <li>• Foot Positioners</li> <li>• Headrest Extension</li> <li>• Tool-less Upholstery Removal</li> <li>• Zippered Pockets in Seat Bottom Upholstery</li> <li>• Elevating Legrest</li> <li>• Made in the USA</li> </ul> |
|---|--|

### Standard Models:

	<u>Price</u>
<input type="checkbox"/> Rodeo Transit RD10T - 904860	<b>\$3,307</b>
<input type="checkbox"/> Rodeo Transit RD12T - 904861	<b>\$3,307</b>
<input type="checkbox"/> Rodeo Transit RD14T - 904862	<b>\$3,307</b>
<input type="checkbox"/> Rodeo Transit RD16T - 904863	<b>\$3,307</b>

### Upholstery Colors

	<u>Price</u>
<input type="checkbox"/> Royal Blue	<b>\$0</b>
<input type="checkbox"/> Forest Green	<b>\$0</b>
<input type="checkbox"/> Cherry Red	<b>\$0</b>
<input type="checkbox"/> Sassy Purple	<b>\$0</b>
<input type="checkbox"/> Princess Pink	<b>\$0</b>
<input type="checkbox"/> Panther Black	<b>\$0</b>
<input type="checkbox"/> Chocolate Brown	<b>\$0</b>
<input type="checkbox"/> Aqua Blue	<b>\$0</b>
<input type="checkbox"/> Apple Green	<b>\$0</b>

RD10	RD12	RD14	RD16
903850	903850	903850	903850
903857	903857	903857	903857
903855	903855	903855	903855
903851	903851	903851	903851
903759	903759	903759	903759
903854	903854	903854	903854
904189	904189	904189	904189
904753	904753	904753	904753
904754	904754	904754	904754

### Transit Required Items

	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Transit Option (Manufacturer Installation Required)	<b>\$0</b>	
<input type="checkbox"/> Headrest Extension	<b>\$0</b>	<b>E0966</b>
<input type="checkbox"/> Headrest Extension Hardware (Required for Headrest Extension)	<b>N/C</b>	
<input type="checkbox"/> H-Harness with Padded Covers	<b>\$0</b>	<b>E0960</b>
<input type="checkbox"/> 3-Pt Positioning Belt	<b>\$0</b>	<b>E0978</b>
Foot Positioners - Select Size Below:		
<input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair)	<b>\$0</b>	<b>K0108</b>
<input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair)	<b>\$0</b>	<b>K0108</b>
<input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair)	<b>\$0</b>	<b>K0108</b>

RD10	RD12	RD14	RD16
904341B	903167B	903168B	903169B
903332B	902863B	900484B	900138B
905462	905462	905462	905462
903988B	903988B	903988B	903988B
904017B	904017B	904017B	904017B
903415B	903415B	903415B	903415B
903416B	903416B	903416B	903416B
903417B	903417B	903417B	903417B

### Transit Options

	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Q'Straint Transit Lap Belt	<b>\$92</b>	<b>E0978</b>
<input type="checkbox"/> Sure-Lok Transit Lap Belt	<b>\$92</b>	<b>E0978</b>
<input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Small ^	<b>\$155</b>	<b>E0950</b>
<input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Medium ^	<b>\$155</b>	<b>E0950</b>

^ Crash tested Foam Tray with vinyl cover. Flip-Up Armrests required.

RD10	RD12	RD14	RD16
903094	903094	903094	903094
903232	903232	903232	903232
903669	903669		
		903670	903670

### Pelvic Positioning

	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Support - General Use Cushion * <sup>1</sup>	<b>\$247</b>	<b>E2601</b>
<input type="checkbox"/> Position Cushion - Medial Thigh Support/Anti Thrust * <sup>1</sup>	<b>\$478</b>	<b>E2607</b>
<input type="checkbox"/> Align Cushion - Lateral Pelvis/Lateral Thigh Support * <sup>1</sup>	<b>\$478</b>	<b>E2607</b>
<input type="checkbox"/> Incontinence Cushion Liner (Cushion req) <sup>1</sup>	<b>\$74</b>	<b>E2619</b>

\* Solid Seat Insert included with all Cushions.  
<sup>1</sup> Manufacturer Installation Required.

RD10	RD12	RD14	RD16
904347	904350	904353	904356
904348	904351	904354	904357
904349	904352	904355	904358
902324	902324	902324	902324

### Foot & Leg Positioning

	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Individual Angle Adjustable Footplates <sup>1</sup>	<b>\$183</b>	<b>K0040</b>
<input type="checkbox"/> Full Padded Footbox (One-Piece Footplate Required)	<b>\$333</b>	<b>K0108</b>
<input type="checkbox"/> Legrest Extension - 2" * <sup>1</sup>	<b>\$56</b>	
<input type="checkbox"/> Medial Thigh Support (Abductor)	<b>\$166</b>	<b>E0957</b>
<input type="checkbox"/> Lateral Thigh Support (Adductor)	<b>\$166</b>	<b>E0956</b>
<input type="checkbox"/> Calf Panel	<b>\$63</b>	<b>E0995</b>
<input type="checkbox"/> Padded Elevating Legrest Cover	<b>\$15</b>	

\* 2" Legrest Extension interferes with caster swivel when seat is in upright position.  
Seat must be tilted back at least 5 degrees OR Legrests must be elevated to allow caster swivel.  
<sup>1</sup> Manufacturer installation required if model was made prior to June 2019.

RD10	RD12	RD14	RD16
	904957	904958	904959
904953	904954	904955	904956
904095	904095	904096	904097
904359	902667	902647	903140
904360	902668	902648	903141
904361	904361	904361	904361
904619	904619	904619	904619

<u>Trunk Positioning</u>	<u>Price</u>	<u>HCPCS</u>
Right Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req)		
<input type="checkbox"/> X Small= 2.75"x 3.75" *	\$121	E0956
<input type="checkbox"/> Small = 3"x 5"	\$121	E0956
<input type="checkbox"/> Medium = 3.75"x 5.5"	\$121	E0956
<input type="checkbox"/> Large = 4.75"x 6"	\$121	E0956
Left Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req)		
<input type="checkbox"/> X Small= 2.75"x 3.75" *	\$121	E0956
<input type="checkbox"/> Small = 3"x 5"	\$121	E0956
<input type="checkbox"/> Medium = 3.75"x 5.5"	\$121	E0956
<input type="checkbox"/> Large = 4.75"x 6"	\$121	E0956
<input type="checkbox"/> Right S/A Adj. Hardware for Lateral Trunk Support	\$244	E1028
<input type="checkbox"/> Left S/A Adj. Hardware for Lateral Trunk Support	\$244	E1028
<input type="checkbox"/> Full Torso Support Vest	\$244	
<input type="checkbox"/> Reducer Seat Insert (Reduces Seat Width 2" & Seat Depth 1")	\$362	

\* RD10 requires X Small laterals.

<u>Head Supports</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Padded Headwings	\$250	E0955
<input type="checkbox"/> Occi Headwings	\$234	E0955

<u>Frame Option</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Recline Option (Requires Manufacturer Installation)	\$473	E1225
<input type="checkbox"/> Height Adjustable Flip-Up Armrests (pair)	\$234	E0973
<input type="checkbox"/> Solid Side Panels	\$216	
<input type="checkbox"/> Mesh Side Panels	\$152	

<u>Tire Options</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> 12.5" Rear Solid Knobby Tire	\$204	
<input type="checkbox"/> 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tire	\$306	

<u>Accessories</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Headrest Cover (Canopy)	\$234	
<input type="checkbox"/> Extended Headrest Cover - No Windows (Canopy)	\$256	
<input type="checkbox"/> Extended Headrest Cover - with Windows (Canopy)	\$302	
<input type="checkbox"/> Upper Extremity Support Surface (Clear Tray) ^	\$114	E0950
<input type="checkbox"/> Adj. Removable Hardware (Required for Clear Tray)	\$233	E1028
<input type="checkbox"/> Medical Necessity Storage Basket	\$302	
<input type="checkbox"/> Under Seat Storage Basket - Tilt-In-Space	\$77	
<input type="checkbox"/> IV Pole (height adjustable and collapsible) <sup>2</sup>	\$154	K0105
<input type="checkbox"/> Oxygen Tank Bag	\$183	E2208
<input type="checkbox"/> Utility Bag	\$39	
<input type="checkbox"/> Attendant Hand Brakes (pair) <sup>1</sup>	\$242	
<input type="checkbox"/> Curb Tipper	\$100	
<input type="checkbox"/> Caster Lock for 2" Wide Front Tires <sup>1</sup>	\$96	K0073
<input type="checkbox"/> Rear Anti-Tip Tubes (Pair)*	\$118	E0971
<input type="checkbox"/> Rain Cover Size 1 (Canopy Req)	\$154	
<input type="checkbox"/> Rain Cover Size 2 (Canopy Req) (RD14 Extended Canopy Req)	\$154	
<input type="checkbox"/> Rain Cover Size 3 (Extended Canopy Req)	\$154	
<input type="checkbox"/> Mosquito Net Size 1 (Canopy req)	\$154	
<input type="checkbox"/> Mosquito Net Size 2 (Canopy req) (RD14 Extended Canopy Req)	\$154	
<input type="checkbox"/> Mosquito Net Size 3 (Extended Canopy req)	\$154	

^ Flip-Up Armrests and Hardware required.  
 \*Highly recommended for users weighing 50lbs or less.  
<sup>1</sup> Manufacturer Installation Required.  
<sup>2</sup> Anti-tip Tubes Required

RD10	RD12	RD14	RD16
904385	904385	904385	904385
	904071	904071	904071
	904072	904072	904072
	904073	904073	904073
904465	904465	904465	904465
	904466	904466	904466
	904467	904467	904467
	904468	904468	904468
904362	904362	904362	904362
904363	904363	904363	904363
904364	903261	903260	903259
904211	904211	904110	904110

RD10	RD12	RD14	RD16
904365	902673	902653	903147
905270	905270	905271	905271

RD10	RD12	RD14	RD16
904342	904343	904344	904345
904346	904346	904346	904346
904388	902672	902652	903145
904389	902680	902681	903146

RD10	RD12	RD14	RD16
904117	904117	904117	904117
904118	904118	904118	904118

RD10	RD12	RD14	RD16
904366	902679	902659	903152
904728	904729	904730	904731
904744	904745	904746	904747
904367	903405	903401	903402
903369	903369	903369	903369
904476	904476	904477	904477
904529	904529	904530	904530
904566	904566	904566	904566
902870	902870	902870	902870
904531	904531	904531	904531
		904370	904371
905310	905310	905310	905310
904650	904650	904650	904650
904649	904649	904649	904649
904939	904939	904939	
		904940	904940
			905283
904655	904655	904655	
		904656	904656
			905319

The Rodeo Transit option includes four manufacturer installed red wheelchair transportation anchors. The Rodeo Transit model has been crash tested and conforms to RESNA WC19 standards. Anthropomorphic Test Device weights are RD10T, RD12T 66lbs/30kg, RD14T 100lbs/45kg, and RD16T 140lbs/64kg.

The Convaid Rodeo mobility base was tested with Rodeo seating module equipped with head support, H-Harness with padded Covers, 3-point restraint system, foot positioners, and 4 transit anchor brackets.

It is Convaid's position that the wheelchair user should be transferred from the wheelchair into appropriate vehicle seating whenever possible. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair, therefore Convaid cannot recommend a transit system.

The WC19 standard is specially designed for the wheelchair to be secured in a forward facing configuration and not in any other direction. For more information on RESNA WC-19 and WC-20 standards please visit <http://www.rercwts.pitt.edu/WC19.html>.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

**For expedited shipping rates, please call 1-844-876-6245**  
 2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166  
 Prices in USD & subject to change, ex works Torrance, California.

Convaid  
 2830 California Street  
 Torrance, CA 90503  
 1 844-US-MOBILITY (844-876-6245)  
 www.etac.us.com

Effective: 5/1/2021

OF0108-ENG-REV.05

