

Convoid Rodeo Tilt w/out Seating - HCPCS E1234

Retail Price List / Order Form

Account #: _____

Purchase Order #: _____

Date: _____

Contact: _____

Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

Rodeo Tilt Standard Features:

- 5° to 45° of Adjustable Tilt
- Silver Powder Coated Frame
- Adjustable One-Piece Push Handle
- Elevating Legrest / One-Piece, Angle Adjustable Footplate
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tires
- Quick Release Rolko Rear Wheels
- Single Foot Activated Wheel Lock Bar
- Seat Depth Growth of 5"
- Made in the USA

Standard Models:

	<u>Price</u>
<input type="checkbox"/> Rodeo w/out Seating RD10WOS - 904963	\$2,376
<input type="checkbox"/> Rodeo w/out Seating RD12WOS - 904964	\$2,376
<input type="checkbox"/> Rodeo w/out Seating RD14WOS - 904965	\$2,376
<input type="checkbox"/> Rodeo w/out Seating RD16WOS - 904966	\$2,376

- Seating Module
 Depth Adjustable Solid Seat Pan (Manufacturer Installation Required)

Price HCPCS
\$154 **E2293**

RD10	RD12	RD14	RD16
904971	904972	904973	904974

- Transit Required Items
 Transit Option (Manufacturer Installation Required)
 Using other manufacturer's Head Support
 2-Point Positioning Belt

Price HCPCS
\$205
\$0
\$46 **E0978**

RD10	RD12	RD14	RD16
904341	903167	903168	903169
995681	995681	995681	995681
903931	903931	904559	904559

- Transit Options
 Q'Straint Transit Lap Belt
 Sure-Lok Transit Lap Belt
 Upper Extremity Support Surface (Foam Tray) Small ^
 Upper Extremity Support Surface (Foam Tray) Medium ^
 ^ Crash tested Foam Tray with vinyl cover. Flip-Up Armrests required.

Price HCPCS
\$92 **E0978**
\$92 **E0978**
\$155 **E0950**
\$155 **E0950**

RD10	RD12	RD14	RD16
903094	903094	903094	903094
903232	903232	903232	903232
903669	903669		
		903670	903670

- Foot & Leg Positioning
 Individual Angle Adjustable Footplates
 Full Padded Footbox (One-Piece Footplate Required)
 Legrest Extension - 2" *
 Foot Positioners - Small 7" - 10" (pair)
 Foot Positioners - Medium 8" - 11" (pair)
 Foot Positioners - Large 9" - 12" (pair)

Price HCPCS
\$183 **K0040**
\$333 **K0108**
\$56
\$160 **K0108**
\$160 **K0108**
\$160 **K0108**

RD10	RD12	RD14	RD16
	904957	904958	904959
904953	904954	904955	904956
904095	904095	904096	904097
903415	903415	903415	903415
903416	903416	903416	903416
903417	903417	903417	903417

* 2" Legrest Extension interferes with caster swivel when seat is in upright position.
 Seat must be tilted back at least 5 degrees OR Legrests must be elevated to allow caster swivel.

- Frame Option
 Recline Option (Manufacturer Installation Required)
 Height Adjustable Flip-Up Armrests (pair)
 Solid Side Panels
 Mesh Side Panels

Price HCPCS
\$473 **E1225**
\$234 **E0973**
\$216 **K0108**
\$152 **K0108**

RD10	RD12	RD14	RD16
904342	904343	904344	904345
904346	904346	904346	904346
904388	902672	902652	903145
904389	902680	902681	903146

- Tire Options
 12.5" Rear Solid Knobby Tire
 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tire

Price HCPCS
\$204
\$306

RD10	RD12	RD14	RD16
904117	904117	904117	904117
904118	904118	904118	904118

- Canopies
 Headrest Cover (Canopy)
 Extended Headrest Cover - No Windows (Canopy)
 Extended Headrest Cover - with Windows (Canopy)
 Canopy Upholstery (Please Choose One Below)
 Royal Blue
 Forest Green
 Cherry Red
 Sassy Purple
 Princess Pink
 Panther Black
 Chocolate Brown
 Aqua Blue
 Apple Green

Price HCPCS
\$234
\$256
\$302
\$0
\$0
\$0
\$0
\$0
\$0
\$0
\$0
\$0
\$0
\$0

RD10	RD12	RD14	RD16
904366	902679	902659	903152
904728	904729	904730	904731
904744	904745	904746	904747
903850	903850	903850	903850
903857	903857	903857	903857
903855	903855	903855	903855
903851	903851	903851	903851
903759	903759	903759	903759
903854	903854	903854	903854
904189	904189	904189	904189
904753	904753	904753	904753
904754	904754	904754	904754

Accessories	Price	HCPCS
<input type="checkbox"/> Upper Extremity Support Surface (Clear Tray) ^	\$114	E0950
<input type="checkbox"/> Adj. Removable Hardware (Required for Clear Tray)	\$233	E1028
<input type="checkbox"/> Medical Necessity Storage Basket	\$302	
<input type="checkbox"/> Under Seat Storage Basket - Tilt-in-Space	\$77	
<input type="checkbox"/> IV Pole (height adjustable and collapsible) ²	\$154	K0105
<input type="checkbox"/> Oxygen Tank Bag	\$183	E2208
<input type="checkbox"/> Utility Bag	\$39	
<input type="checkbox"/> Attendant Hand Brakes (pair) ¹	\$242	
<input type="checkbox"/> Curb Tipper	\$100	
<input type="checkbox"/> Caster Lock for 2" Wide Front Tires ¹	\$96	K0073
<input type="checkbox"/> Rear Anti-Tip Tubes (Pair)*	\$118	E0971
<input type="checkbox"/> Rain Cover Size 1 (Canopy Req)	\$154	
<input type="checkbox"/> Rain Cover Size 2 (Canopy Req) (RD14 Extended Canopy Req)	\$154	
<input type="checkbox"/> Rain Cover Size 3 (Extended Canopy Req)	\$154	
<input type="checkbox"/> Mosquito Net Size 1 (Canopy req)	\$154	
<input type="checkbox"/> Mosquito Net Size 2 (Canopy req) (RD14 Extended Canopy Req)	\$154	
<input type="checkbox"/> Mosquito Net Size 3 (Extended Canopy req)	\$154	

RD10	RD12	RD14	RD16
904367	903405	903401	903402
903369	903369	903369	903369
904476	904476	904477	904477
904529	904529	904530	904530
904566	904566	904566	904566
902870	902870	902870	902870
904531	904531	904531	904531
		904370	904371
905310	905310	905310	905310
904650	904650	904650	904650
904649	904649	904649	904649
904939	904939	904939	
		904940	904940
			905283
904655	904655	904655	
		904656	904656
			905319

^ Flip-Up Armrests and Hardware required.
 *Highly recommended for users weighing 50lbs or less.
¹ Manufacturer Installation Required.
² Anti-tip Tubes Required.

The Rodeo Transit option includes four manufacturer installed red wheelchair transportation anchors. The Rodeo Transit model has been crash tested and conforms to RESNA WC19 standards. Anthropomorphic Test Device weights are RD10T, RD12T 66lbs/30kg, RD14T 100lbs/45kg, and RD16T 140lbs/64kg.

The Convaid Rodeo Mobility Base was tested with the Rodeo Seating Module, equipped with a head support and wheelchair-anchored 2-point pelvic lap belt. Do not use the foot positioners on the Convaid Rodeo during transit in a motor vehicle. The performance and safety of accessories and/or other seating options not tested by Convaid is unknown, and no warranty or claim about their performance is or can be made.

It is Convaid's position that the wheelchair user should be transferred from the wheelchair into appropriate vehicle seating whenever possible. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair, therefore Convaid cannot recommend a transit system.

The WC19 standard is specifically designed for the wheelchair to be secured in a forward facing configuration and not in any other direction. For more information on RESNA WC-19 and WC-20 standards please visit <http://www.rercwts.pitt.edu/WC19.html>.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245
 Prices in USD & subject to change, ex works Torrance, California.

Effective: 5/1/2021

Convaid
 2830 California Street
 Torrance, CA 90503
 1 844-US-MOBILITY (844-876-6245)
 www.etac.us.com

